




FUJIFILM

Value from Innovation

LAL Reagent Products for Detection of BACTERIAL ENDOTOXIN



Wako

Greetings from FUJIFILM Wako Chemicals U.S.A. Corp.

Since the establishment of our first satellite sales office in Dallas, TX (1981), to the construction of our corporate headquarters and manufacturing facilities in Richmond, VA (1989), Fujifilm Wako Chemicals U.S.A. Corporation has strived to provide customers in all scientific disciplines with products of the utmost quality and dependability.

Long recognized as a world-renowned supplier of high purity chemicals and reagents, our company continues to maintain a proud history of product quality and customer service through the establishment of the LAL Division, and the introduction of our new PYROSTAR™ ES-F line for the detection of bacterial endotoxin.

This publication represents the culmination of more than 30 years of research and development, dedicated to providing our customers with endotoxin-specific reagents for “every user and for every method”.

We invite you to review our catalog and look forward to having the opportunity to serve you.

Our Promise

As an FDA licensed facility, FUJIFILM Wako Chemicals U.S.A. Corp. – LAL Division is committed to ensuring that our production site and LAL reagents comply with all the rules, regulations, and quality standards set forth by FDA for current Good Manufacturing Practices (cGMP's).

Horseshoe Crab Conservation

FUJIFILM Wako Chemicals U.S.A. Corp. is very much concerned with maintaining the viability of the horseshoe crab population. We are dedicated to following practices that ensure the careful handling and good quality of crabs used for LAL manufacture that both minimize injury and protect this invaluable species. After bleeding, the crabs are returned the next day by our fisherman to the same waters where they were collected. To assist in the collection of data for crab conservation studies, FUJIFILM Wako Chemicals U.S.A. Corp. participates in a horseshoe crab tagging and monitoring program coordinated by the U.S. Fish and Wildlife Service.



FUJIFILM Wako Chemicals U.S.A. Corp. - Cape Charles Bleeding Facility: Cape Charles, Virginia 23310

Endotoxin-Specific LAL Reagents	2
PYROSTAR™ ES-F Series	2-4
Limulus Color KY Series	5
Limulus PS Single Test	6
Toxinometer® Measurement System	7
Endotoxin Test Related Accessories	8-11
FUJIFILM Wako Distributors	12-13

Wako Worldwide Locations

United States of America

FUJIFILM Wako Chemicals U.S.A. Corp.
1600 Bellwood Road
North Chesterfield, VA 23237
Toll Free: 800-992-9256
Tel: 804-714-1919
Tel: 804-672-4655
Fax: 804-271-7791
Email: wkuspyrostarinfo@fujifilm.com
www.wakopyrostar.com

Japan

FUJIFILM Wako Pure Chemical Corp.
1-2 Doshomachi
3-Chome, Chuo-Ku,
Osaka 540-8605
Japan
Telephone: + 81-6-6203-3741
Facsimile: + 81-6-6222-1203
www.wako-chem.co.jp

Germany

FUJIFILM Wako Chemicals Europe GmbH
Fuggerstraße 12
41468 Neuss
Germany
Telephone: + 49-2131-311-0
Facsimile: + 49-2131-311-100
www.wako-chemicals.de



FUJIFILM Wako Chemicals U.S.A. Corp introduces Endotoxin-Specific LAL Reagents for every user and for every method

Untreated LAL reagents react not only with endotoxin but also with (1→3) β-D-glucan, a fungal cell wall component that initiates the clotting cascade by activating the Factor G pathway. To eliminate this (1→3) β-D-glucan activation, various endotoxin-specific reagents have been developed by either removing Factor G or inhibiting its activation.

FUJIFILM Wako Chemicals U.S.A. Corp now provides a new line of reagents formulated to perform each of the three endotoxin detection assays, namely the Gel-Clot Assay, Kinetic-Turbidimetric Assays (KTA) and Kinetic-Chromogenic Assays (KCA).

The PYROSTAR™ ES-F Series

The PYROSTAR™ ES-F series of reagents represent our uniquely endotoxin specific line of Limulus Amebocyte Lysate products, which are formulated to be used either qualitatively as a gel-clot reagent or quantitatively as a Kinetic-Turbidimetric reagent.

Product features

- Endotoxin-specific lysate, avoids false positive results from glucans
- Available in multi-tests vials or single-test vials
- Can be used as either a gel-clot or Kinetic-Turbidimetric Assay (KTA) reagent
- (KTA) assays can be performed in tube reader or microplate reader
- Gel-Clot lysate sensitivities range from 0.015 to 0.25 EU/mL
- Available with matched control standard endotoxin (CSE)
- PYROSTAR™ ES-F reagents are available with a KTA quantitative range of either 0.001 EU/mL to 10 EU/mL.
The KTA quantitative range is related to the Gel-Clot sensitivity.
- 100ul sample size when used with tube reader; 50ul sample size when used with microplate reader



Gel-Clot Sensitivity (EU/mL)	KTA Quantitative Range (EU/mL)
0.03 - 0.25	0.01 to 10
0.015	0.001 to 10

Multi-Test Lysate with matched CSE

PYROSTAR™ ES-F 2.0 mL Test Kit • 4 multi-test vials (2.0 mL) + 1 vial CSE (500 ng/vial)		
Catalog number	Gel-clot Sensitivity EU/mL	Kinetic-Turbidimetric Assay Quantitative Range (EU/mL)
WPEK4-20015	0.015	0.001 to 10
WPEK4-20003	0.03	0.01 to 10
WPEK4-20006	0.06	
WPEK4-20125	0.125	
WPEK4-20025	0.25	



PYROSTAR™ ES-F 5.2 mL Test Kit • 4 multi-test vials (5.2 mL) + 1 vial CSE (500 ng/vial)		
Catalog number	Gel-clot Sensitivity EU/mL	Kinetic-Turbidimetric Assay Quantitative Range (EU/mL)
WPEK4-50015	0.015	0.001 to 10
WPEK4-50003	0.03	0.01 to 10
WPEK4-50006	0.06	
WPEK4-50125	0.125	
WPEK4-50025	0.25	



Single- Test Lysate with matched CSE

PYROSTAR™ ES-F Single Test Kit • 25 single-test vials + 1 vial CSE (500 ng/vial)		
Catalog number	Gel-clot Sensitivity EU/mL	Kinetic-Turbidimetric Assay Quantitative Range (EU/mL)
WPESK-0015	0.015	0.001 to 10



Bulk Lysate

PYROSTAR™ ES-F 2.0 ml Bulk Kit • 100 multi-test vials (2.0 mL)		
Catalog number	Gel-clot Sensitivity EU/mL	Kinetic-Turbidimetric Assay Quantitative Range (EU/mL)
WPEM-20015	0.015	0.001 to 10
WPEM-20003	0.03	0.01 to 10
WPEM-20006	0.06	
WPEM-20125	0.125	
WPEM-20025	0.25	



Bulk Lysate

PYROSTAR™ ES-F 5.2 ml Bulk Kit • 100 multi-test vials (5.2 mL)		
Catalog number	Gel-clot Sensitivity EU/ml	Kinetic-Turbidimetric Assay Quantitative Range (EU/mL)
WPEM-50015	0.015	0.001 to 10
WPEM-50003	0.03	0.01 to 10
WPEM-50006	0.06	
WPEM-50125	0.125	
WPEM-50025	0.25	



The Limulus Color KY Series

The Limulus Color KY Series includes both a multi-test kit and a single-test kit, each designed for time-based chromogenic analysis, using a synthetic substrate which produces a yellow color and can specifically detect endotoxin with high sensitivity.



Product features

- Endotoxin-specific lysate, avoids false positive results from glucans
- Available in multi-tests vials or single-test vials
- Quantitative Kinetic-Chromogenic Assay (KCA) reagent
- KCA assays can be performed in tube reader or microplate reader
- KCA quantitative range detection limit of 0.0002 EU/mL (single-type) and 0.0005 EU/mL (multi-type).
- Available with matched control standard endotoxin (CSE)
- 100ul sample size when used with tube reader; 50ul sample size when used with microplate reader

Lysate with matched CSE

Limulus Color KY Test Kit • 3 multi-test vials (2.0 mL) + 1 vial CSE (500 ng/vial)	
Catalog number	Number of Tests
291-53101	60 tests

Limulus Color KY Single Test Kit • 25 single test vials + 1 vial CSE (500 ng/vial)	
Catalog number	Number of Tests
291-53601	25 tests

Limulus PS Single Test

The Limulus PS Single Test is a kit composed of LAL endotoxin specific reagent and an affinity resin suspension known as PyroSep™, designed to overcome any product interference by adsorbing potential endotoxin in samples while washing away the inhibitory components.



Product Features

- Endotoxin-specific lysate, avoids false positive results from glucans
- KTA assay is performed in tube reader
- Allows endotoxin specific measurement of samples which inhibit traditional LAL assays
- Fat-soluble vitamins, oils and fats can be tested if they can be dissolved in ethanol.

Limulus PS Single Test	
Catalog Number	Contents
299-54501	20 LAL ES single test vials PyroSep™ resin suspension LAL reconstitution solution Wash solution Sample diluent 20 Glass capillary columns

PS Accessory Kit (optional)	
Catalog Number	Contents
294-33311	12 syringes (20 mL) 12 T-shape stopcocks R Type 20 Dedicated adapters 12 Bulldog clips 25 Polystyrene tubes (5 mL) 25 Polystyrene tubes (14 mL) Test tube rack

Toxinometer® Measurement System **FUJIFILM**

Value from Innovation

Next generation ET-7000 to be released 4th Quarter of 2018



Toxinometer® ET-6000 - Endotoxin Measurement System

The Toxinometer® ET-6000 is our new computer-operated kinetic incubating tube reader, which is exceptionally user-friendly and easily expandable. Depending on the number of samples to be processed, our state of the art expansion modules can be connected to allow for endotoxin testing in a wide range of fields and sample quantities.

Product Features

- Kinetic Incubating Tube Reader
- Single-test configuration avoids “Hot Wells”
- For use in Kinetic-Turbidimetric, Kinetic-Chromogenic, and Gel-Clot assays
- A single assay module can simultaneously measure up to 16 samples
- Expansion modules are available to extend number of samples in multiples of 16
- Temperature settings at both 30°C and 37°C
- Endotoxin determination in compliance with FDA guidelines as well as with Pharmacopeal (USP/EP/JP) BET monographs
- Toximaster QC7 - State of the art software
 - FDA 21 Part 11 Compliant
 - Allows trending of product results over time
 - Allows early detection of potential product failures
 - Generates hard copy printouts documenting assay parameters and sample results
 - Test to completion without intervention
 - Develop automated assay systems with ease
 - Offers improved precision and accuracy over the gel-clot method



Catalog No.	Model	Power Source	Contents
293-33509	Toxinometer® ET-6000/U Part 11 Set	100-120 +/- 10% VAC (USA)	1 Toxinometer® ET-6000 Toximaster QC7 Software w/5 User Licenses
290-33519	Toxinometer® ET-6000/E Part 11 Set	220-240 +/- 10% VAC (Europe)	1 Personal Computer System Validation Documents
293-33989	Toxinometer® ET-6000/U Non- Part 11 Set	100-120 +/- 10% VAC (USA)	1 Toxinometer® ET-6000 Toximaster QC7 Software w/1 User License
299-33969	Toxinometer® ET-6000/E Non- Part 11 Set	220-240 +/- 10% VAC (Europe)	1 Personal Computer
297-33529	Toxinometer® ET-6000/U Expansion Module	100-120 +/- 10% VAC (USA)	16 well expansion unit
294-33539	Toxinometer® ET-6000/E Expansion Module	220-240 +/- 10% VAC (Europe)	16 well expansion unit

Control Standard Endotoxin

Our Control Standard Endotoxin is a lyophilized product composed of 500ng/vial of endotoxin, purified from E. Coli UKT-B strain that serves as an accurate and dependable standard for endotoxin testing.



Product Features

- Endotoxin derived from E.Coli UKT-B
- Can be used to prepare controls and standard curves
- RSE/CSE ratios supplied to match specific lots of LAL
- Reconstituted CSE can be stored at 2-10°C for 1 month

(Ready-made indicators and new high concentration LPS to launch Fall of 2018)

Catalog Number	Contents
CSE4037-5006	6 vials (500 ng/vial)



Endotoxin Test Related Accessories



BioClean® Series

FUJIFILM Wako's premiere collection of endotoxin-free pipette tips, test tubes and caps designed for worry free endotoxin testing on our Toxinometer® measurement system.

Product Features

- Endotoxin-free pipette tip (<0.005 EU/tip)
- Individually packaged, ideal for use in clean rooms
- Packaging film features dust-free protection
- Sterilized by gamma irradiation (25kg ray)
- Tips fit Finnpiquette digital micro pipetters (20-200ul or 100-1000ul sizes)
- Endotoxin-free tubes and caps sterilized at 250 °C
- For use with Toxinometer® measurement system



Endotoxin-free Pipette Tips and Microplates

Catalog Number	Product Name	Volume Size	Quantity
291-35021	BioClean® Tip Wako 200 II	200 µl	100 tips
298-35031	BioClean® Tip Wako 1000 II	1000 µl	100 tips
294-35011	BioClean® Tip Wako Extend S II	200 µl	100 tips
293-35221	Bioclean® Microplate	96 well	50 pcs/pk



Endotoxin-free Test Tubes and Caps

Catalog Number	Product Name	Dimensions	Quantity
292-32751	Limulus Test Tube-S w/Aluminum Cap	12 x 75mm	80 pcs
293-26551	Limulus Test Tube-S	12 x 75mm	100 pcs
293-28251	Aluminum Cap-S	14.7 x 18mm	100 pcs

Gel-clot Reaction Tubes

Our Gel-clot Reaction Tubes are made of high quality Borosilicate glass and certified free of interfering endotoxins and glucans.

Product Features

- Endotoxin-free test tubes (<0.001 EU/tube)
- Designed for Gel-clot testing
- Made with Borosilicate glass



Catalog Number	Product Name	Quantity
CT-1075	Gel-clot Reaction Tube	200 tubes/pk
CT-1075C	Gel-clot Reaction Tube w/ Aluminum Cap	50 pcs/pk
CT-1075B	Gel-clot Reaction tube, Bulk	1250 tubes/pk
DL-13100	Depyrogenated Dilution Tube 13 x 100mm	50 pcs/pk

LAL Reagent Water

Wako proudly provides high-quality, endotoxin free water ideal for all of your LAL testing needs.

Product Features

- Endotoxin-free water (<0.001 EU/mL)
- Steam sterilized by USP standards
- Derived from water for injection, USP
- Non LAL reactive

(125mL in plastic bottles to launch Fall 2018)



Catalog Number	Product Name	Quantity
LRW-12100	Lysate Reagent Water, 100ml	12 x 100mL vials
LRW-2030	Lysate Reagent Water, 30ml	20 x 30mL vials

Endotoxin Extracting Solution

FUJIFILM
Value from Innovation

Endotoxin Extracting Solution for LAL Test

Traditionally, water or saline solution has been used to extract endotoxin in tests on medical devices and equipment; however, the efficacy of this extraction method has recently come into question. In order to provide our customers with a more reliable method for medical device testing, FUJIFILM Wako has developed an endotoxin-free extracting solution containing human serum albumin (HSA) that is capable of extracting endotoxins which cannot be extracted in water or saline solution.

Product Features

- Capable of extracting endotoxins from surfaces which cannot be extracted in water or saline
- Recommended for use in endotoxin testing on equipment and devices which may come into contact with fluids containing blood or protein.



Catalog Number	Contents
293-51601	4 x 10mL vials



Contact Our Distributors

Looking for a distributor in your region?
Please see our list of distributors around the globe:

DOMA INTERNATIONAL LLC

(Distributor for Argentina, Paraguay, Ecuador, Chile, and Bolivia)

5717 NW 125 Ave.

Coral Springs

Florida, 33076

USA

Contact Person: Lic. Marcelo Madeo

Email: info@domainternationalllc.com

ENDOLAB S.A. (Distributor for Argentina)

Cuba 3812, Piso 8 Oficina "B",

Capital Federal, Buenos Aires,

Argentina

Contact Person: Mr. Roberto Aunon Elliff

Phone: (54-11) 4701-7734

Email: gerencia@endolabsa.com.ar

Website: www.endolabsa.com.ar

Lobov Científica Div. of Lobov y Cía. S.A.C.I.

(Distributor for Argentina)

Av. Franklin D. Roosevelt 5828

CABA, C1431BZV

Argentina

Phone: +54-11 4574-3500 / Fax: 4574-3600

Fax: +54-11 4574-3600

Email: lobov@lobov.com.ar

Website: www.lobov.com.ar

Bioaid (Distributor for Bangladesh)

House# 1/B/1, Unit-A1, 1st Floor, Road#01

Shyamoli, Mohammadpur

Dhaka-1207

Bangladesh

Phone: +88029134673

Phone: 01965864733

M/S RIFA Enterprise (Distributor for Bangladesh)

51 North Goran, Dhaka-1219

Bangladesh

Sales Office:

1 New Eskaton Road,

Razzak Plaza, 10th Floor, Dhaka-1000, Bangladesh

Phone: +8801911-367920

Phone: +880-2- 9342501

DP Unilab (Distributor for Belarus and Russia)

220094, Minsk, Gorovtza str.24, app.3

Republic of Belarus

Contact Person: Eugene Kostoglod

Phone: +375 17 3195791

Fax: +375 17 3195791

Email: dpunilab@telecom.by

Website: www.unilab.by

Sellex S.A.C. (Distributor for Brazil)

Rua Arandu, 205 / 1105,

São Paulo, SP,

Brazil

Email: vendas@sellex.com

Khema Scientific Co., Ltd. (Distributor for Cambodia)

#5CHE0E1, Road 5CH03, Phum Trapeang Thloeng

Sangkat Chom Chao, Khan Porsenchey

Phnom Penh, Cambodia

Contact Person: Mr. San Sokhem

Phone: +855-(0)17-667311

Email: info@khemasci.com

Cedarlane (Distributor for Canada)

4410 Paletta Court

Burlington, Ontario L7L 5R2

Canada

Contact Person: Kyle Pahapill

Phone: 289-288-0001

Email: sales@cedarlanelabs.com

Website: www.cedarlanelabs.com

Pulse Scientific Inc. (Distributor for Canada)

18 - 5100 South Service Road

Burlington, Ontario, L7L6A5

Canada

Contact Person: Allen Chan

Phone: 905-3338188 or 800-3637907

Email: allen@pulsescientific.com

or order@pulsescientific.com

Website: www.pulsescientific.com

Water Treatment Systems S.A.S

(Distributor for Colombia)

Diagonal 59 # 22a 45, Bogotá

Colombia

Contact Person: Jose Luis Pinzon

Phone: 0571 7322650 /4633 cell 3208802996

Femto Trade Co. (Distributor for Egypt)

Villa No. 25 , Ibn Battuta St. Haram

Giza, 12111

Egypt

Contact Person: Dr. Essam Abdallah

Phone : 0020233928140

Phone: 00201005126629

Fax: 0020233928424

Email sales@femtotrade.com

Email: Essamfemto@yahoo.com

BINOCULAR Trading (Distributor for Ethiopia)

Bole Sub City Addis Ababa, 19996

Ethiopia

Contact Person: Zenaw Workneh

Phone: +251-911 733 801

African Bio-expertise Services Ltd

(Distributor for Ghana, Ivory Coast, Togo, and Benin)

Boulevard de la 8ème Tranche

Cocody Angré

Abidjan - Côte d'Ivoire

Tel: +225 22 42 89 02 / +225 07 35 28 97

Contact Person: Dr Amadou DIALLO

Email: enquiry@abs-ai.com

Website: www.abs-ai.com

Bio Solutions (Distributor for India)

102, Ashok Heights

Dr. S. Radhakrishnan Marg

Andheri (East), Mumbai-400069

Fax: 022-26824494

Phone: 9819973583

Email: mail@biosolutions.co.in

Contact Us Our goal is to provide you with quality customer service. Whether it is assisting with your order, providing technical information, or advising pricing on our products, the staff of Fujifilm Wako Chemicals U.S.A. Corporation is here to help. Please contact us with your inquiries directly or through our website www.wakopyrostar.com

FUJIFILM Wako Distributors

FUJIFILM
Value from Innovation

PAMMVI EXPORTS PVT LTD (Distributor for India)

C 501 Remi Bizzcourt, Off Veera Desai Road Andher West Mumbai, Maharashtra, 400053
India
Phone: 91 22 26731904 / 26731052 / 26732760
Contact Person: Managing Director, MD Viveck Gupta
Email: viveck.guptaa@pammvi.com
Contact Persons: Technical Director, Amit Gupta
Email: amit.gupta@pammvi.com
Website: www.pammvi.com

Advansys Technology for Life Ltd. (Distributor for Israel)

1 Hamada Street; Science Park
Rehovot, 7670301
Israel
Phone: +972-8-9417118
Fax: +972-8-9492565
Website: www.advansys.co.il

Edo-Elemed (Distributor for Lebanon, Jordan and Syria)

St Takla Boulevard, Jdaydeh, Bldg.# 4, Bierut
Lebanon
Phone: +961 76 63 11 82
Email: info@edo-elemed.com, raymond.rahall@edo-elemed.com
Website: www.edo-elemed.com

BIOTEK DOO SKOPJE (Distributor for Macedonia)

Zagrebska 28V 1000 Skopje, Ecully
Macedonia
Contact Person: Maja Andonova
Phone: 00 389 2 3110106
Email: mandonova@biotek.com.mk; sales@biotek.com.mk
Website: www.biotek.com.mk

Bio-Focus Saintifik Sdn Bhd (Distributor for Malaysia)

13A, Block H, Jalan PJU 1A/3,
Taipan 2 Damansara, Ara Damansara, Petaling Jaya
Selangor Darul Ehsan, 47301 Malaysia
Phone: +603 7842 9850
Email: sales@bio-focus.com

Bioprocesos Farmacéuticos Integral S.A. de C.V. (Distributor for Mexico)

Fontana #28 Col. Frac. Lomas Boulevard
Tlalhepantla de Baz, 54020
Mexico
Phone: (01 55) 15091624, 50121830

Mast Comm (Distributor for Morocco)

2 rue Chenguiti, Et 4 N° 6, Belvedere
Casablanca, 20310
Morocco
Contact Person: Abdelaziz Agouzzi
Phone: 0522249755/66
Email: mastcomm1@gmail.com

Sind Trading (Distributor for Pakistan)

3/14 Arkay Square, Shahrah-e-Liaquat
Karachi, 74000
Pakistan
Phone: 0300-8240578
Email: shujatmohsin@gmail.com
Email: sind_trading_co@yahoo.com

Filter Lab Peru (Distributor for Peru)

Av. Tejada 259 oficina 402 Miraflores
Lima, Peru
Contact Person: Rafael Massa
Phone: 4714596
Phone: 955543528
Email: rmassa@filterlab.com.pe
Website: www.filterlab.com.pe

IFLE CHEMICALS S.A.C (Distributor for Peru)

(Industria Farmacéutica Libre de Endotoxinas Químicas S.A.C.)
Av. Alfonso Ugarte N°1454 OF. 302 - Breña
Lima - Perú
Contact Person: Q.F. Stephany Criollo S.
Phone: (511) 495-7968 / (511) 993-488675
Email: ifle.gerencia@gmail.com / ifle.ventas@gmail.com

B/B Medical LLC (Distributor for Russia)

Office 12, Floor 5, 14, Shukhova str.
Moscow, 115162, Russia
Phone: +74955443902
Contact Person: German Yakovenko Email: sales@bbmed.ru
Website: www.bbmed.ru

Bright Glory Technology Inc. (Distributor for Taiwan)

5F-3, No. 162, Sec. 4
Chung Hsiao E. Road
Taipei, Taiwan
Phone: +886-2-27317727
Fax: +886-2-27318221
Email: bright@bge.com.tw
Website: www.bge.com.tw

Mikromed Medikal (Distributor for Turkey)

Paz. lth. lhr. San. ve Tic. Ltd.Şti.
Mehmet Akif Mah, Elalmiş Cad. Asya Sk, No.19 D.3
34774 Ümraniye, İstanbul – TÜRKİYE
Phone: +90 216 668 00 52/42
Email: info@mikromedltd.com

Friedman-Ukraine Ltd. (Distributor for Ukraine)

Pr. Peremogy, 68/1, r.62,
Kiev, 03113,
Ukraine
Contact Person: Vitaliy Melnik
Phone: +38-067-504-20-79
Email: office@friedman.com.ua

Microgen Medical Equipments Trading (Distributor for United Arab Emirates)

New Industrial Area, Al Ain, Abu Dhabi
United Arab Emirates
Contact Person: Mr. Labeeb Abdulla
Phone: +971506956253
Email: sales@microgen-me.com
Website: www.microgen-me.com

KOMPASS (Distributor for Uruguay)

Prudencio Vazquez y Vega 1112 C.P. 11300, Montevideo
Uruguay
Contact Person: Anibal Bruzzone
Phone: +59827104373
Email: info@kompass.com.uy
Website: http://www.kompass.com.uy

Lab Solutions Ltd.

(Distributor for Zimbabwe, Zambia, and South Africa)

6 Cambitzis Building
Kenny Road, Avondale
Harare, Zimbabwe
Contact Person: Hensley Chamboko
Phone: +263 772 806110
Email: info@labsolutions.co.zw
Website: www.labsolutions.co.zw



FUJIFILM

United States of America

FUJIFILM Wako Chemicals U.S.A. Corp.
1600 Bellwood Road
North Chesterfield, VA 23237
Toll Free: 800-992-9256
Tel: 804-714-1919
Tel: 804-672-4655
Fax: 804-271-7791
Email: wkuspyrostarinfo@fujifilm.com
www.wakopyrostar.com

Japan

FUJIFILM Wako Pure Chemical Corp.
1-2 Doshomachi
3-Chome, Chuo-Ku,
Osaka 540-8605
Japan
Telephone: + 81-6-6203-3741
Facsimile: + 81-6-6222-1203
www.wako-chem.co.jp

Germany

FUJIFILM Wako Chemicals
Europe GmbH
Fuggerstraße 12
41468 Neuss
Germany
Telephone: + 49-2131-311-0
Facsimile: + 49-2131-311-100
www.wako-chemicals.de